(90 × 210mm Size)

〈For Research Use Only〉 Code No. 298-65701 (1,000 Tests)

LabAssay™ Glucose (Mutarotase-GOD method)

For Research Use Only Not for Diagnostic Use

〔Introduction〕

α -D-Glucose and β -D-glucose in solution maintain equilibrium in a constant ratio. Glucose oxidase reacts only with β -D-glucose and does not react with α -D-Glucose. Therefore, α -D-Glucose is converted into β -D-glucose using mutarotase.

LabAssay™ Glucose is the reagent kit for assay of glucose based on an enzymatic method with a combination of mutarotase and glucose oxidase. The kit is used for the quantitative determination of glucose in mouse serum, plasma or urine. It is a simultaneous multi-sample assay format using a microplate, but measurements can also be made using a test tube.

〔Kit contents〕

(1)	Buffer (Phosphate buffer, pH 7.1 60mmol/L Phenol 5.3mmol/L	150mL	2vials
(2)	Chromogen Reagent After reconstitution (Mutarotase 0.13units/mL Glucose oxidase (GOD) 9.0units/mL Peroxidase (POD) 0.65units/mL 4-Aminoantipyrine 0.50mmol/mL Ascorbate oxidase (AOD) 2.7units/mL	For 150mL	2vials
(3)	Glucose Standard I (200mg/dL Glucose)	10mL	1vial
(4)	Glucose Standard II (500mg/dL Glucose)	10mL	1vial

〔Materials and apparatuses to be prepared〕

- 96wells microplate (transparent type)
- Micropipette
- Incubator maintaining at 37℃*
- Plate mixer*
- Microplate reader with 505nm wavelength filter
(* if the microplate reader is not equipped.)



〔For test tube method〕

- Test tube
- Pipette
- Incubator maintaining at 37℃
- Spectrophotometer or colorimeter with 505 nm wavelength filter

—1/10—

〔Assay principle〕

When a sample is mixed with the Chromogen Reagent, the α -form of glucose in the sample is converted to β -form by mutarotase. β -D-Glucose is oxidized and yields hydrogen peroxide by glucose oxidase (GOD). In the presence of peroxidase (POD), the formed hydrogen peroxide yields a red pigment by quantitative oxidation condensation with phenol and 4-aminoantipyrine. The glucose concentration is obtained by measuring absorbance of the red pigment.

〔Preparation of reagents to be used〕

① Chromogen reagent :

Dissolve 1 vial of Chromogen Reagent with 150mL of Buffer. After reconstitution, the solution should be stored at 2-10℃ and used within 1 month.

② Standard solution (Microplate method)

Standard solution is prepared by dilution of the provided Standard.

No.	Glucose Standard I	Glucose Standard II	Distilled or deionized water	Sample volume	Concentration
1	50 μ L	—	150 μ L	2 μ L	50mg/dL
2	100 μ L	—	100 μ L	2 μ L	100mg/dL
3	undiluted	—	—	2 μ L	200mg/dL
4	—	150 μ L	100 μ L	2 μ L	300mg/dL
5	—	200 μ L	50 μ L	2 μ L	400mg/dL
6	—	undiluted	—	2 μ L	500mg/dL

〔Procedure〕

(1) Assay in a microplate

Perform the assay in the wells according to the following table scheme.

	Test	Standard	Blank
Chromogen reagent	300 μ L	300 μ L	300 μ L
Sample	2 μ L	2 μ L	— *1
	Mix well and incubate at 37℃ for 5 min. Measure the absorbance at 505nm*2 of the test sample and standard solution with the blank solution as the control.		

*1 The addition of water is omitted because the absorbance difference in the presence of the addition of 2 μ L of water does not exhibit any practical use.

*2 Use 505 nm (Main) and 600 nm (Sub) for two wavelengths.

—2/10—

(2) Assay in a test tube


Perform the assay in a test tube according to the following table scheme.

	Test	Standard	Blank
Chromogen reagent	3.0mL	3.0mL	3.0mL
Sample	20 μ L	20 μ L	— *1
	Mix well and incubate at 37°C for 5 min. Measure the absorbance of the test sample and standard solution with the blank solution as the control. Colorimeter with 505 nm Spectrophotometer : 505 nm*2		

* 1 The addition of water is omitted because the absorbance difference in the presence of the addition of 20 μ L of water does not exhibit any practical use.

* 2 Use 505 nm (Main) and 600 nm (Sub) for two wavelengths.

[Standard curve] [assay in a microplate]

Microplate reader : SAFIRE (TECAN)

[Performance]

(1) Sensitivity [assay in a test tube]

- The absorbance is below 0.03, when measuring purified water as a sample.
- The absorbance is 0.40 to 0.55, when measuring 200 mg/dL glucose as a sample.

(2) Specificity

The glucose concentration is less than $\pm 12\%$, when measuring the known concentration of a control serum or a urine as a sample.

[Usage Notes]

(1) Sample

- Samples should be used immediately after collecting.
- The glucose in whole blood decreases due to glycolysis in erythrocytes when a sample is kept after collecting blood. Therefore, separation of blood cells should be done as quickly as possible.
- Ascorbic acid and bilirubin cause a slightly negative effect on the assay.
- Hemolysis gives a slightly positive effect on the assay.
- Anticoagulants such as heparin, citrate, oxalate and EDTA and sodium fluoride as a glycolytic inhibitor, should not affect the assay when they are used in their usual amounts.

(2) Notes on the assay

- Do not use the reagents described above in any procedures other than those described herein. Performance cannot be guaranteed if the reagents are used in other procedures or for other purposes.
- Operate the measurement apparatuses according to operator's manuals under appropriate conditions. Consult the apparatus manufacturer for details.
- Keep the provided reagents under the indicated conditions before the expiration.
- Do not use reagents, which were frozen in error. Such reagents may give false results.
- After opening the reagents, it is recommended to use them immediately. The opened vials should be capped and kept under the specified conditions.
- Do not use the containers and other materials in the package for any purposes other than those described herein.
- Do not use the reagents at other reaction temperatures and reaction times than described herein.
- The reaction completes in about 4 min at 37°C. When the incubation is continued, the absorbance gradually decreases. Therefore, the incubation should be done within 20 min.
- Glucose can be also measured by a reaction at room temperature (15°C or more) for 15 min.
- Color level of the reaction solution should show very little change within 1 hour at room temperature.
- This kit is for research use only. Not for diagnostic use.

(3) Safety precautions

- If a reagent comes into contact with the mouth, eyes, or skin, immediately wash with a lot of water. Consult a physician if necessary.
- A pipette with a safety pipette filler should be used.
- Be careful not to cut yourself with the aluminum cap when removing it from the vial.

(4) Waste

- The waste should be processed appropriately according to regulations.
Buffer contains 500 mg/L phenol.
- All the devices including reagents and vials that come into contact with specimens should be considered potentially infectious.

Expiration date : 18 months after the manufacture

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