

PERSONAL PRECISION

NEW & IMPROVED

AN ENHANCED SYSTEM THAT IS FASTER AND MORE FLEXIBLE

IMPROVEMENTS INCLUDE

Extended dye set	More flexibility
Improved signal	Higher sensitivity
Faster run times	Less waiting




EXTENDED DYE SET

New LEDs and optical filters have extended the range of fluorophores that the MyGo Mini S is compatible with, which now include LC Green for HRM or Unlabelled Probe Melting analysis, and ROX for TaqMan assays.


IMPROVED SIGNAL QUALITY


A new CMOS sensor, with increased performance, provides 4.3 x higher signal to noise for challenging applications.



FASTER RUNS

Reduced hold times enable results in less than 30 minutes.







NOW IN RED

The MyGo Mini S is available in red.

All trademarks are the property of their respective owners. Design and specifications are subject to change without notice.

Version 001